

Product Catalogue



About IKON Mobility

IKON Mobility by Redgum Brand has been created through the consistent demand from the supporters of the Redgum Brand range to bring an innovative line up of battery powered electric mobility devices to the market.

The Redgum Brand has now been servicing the needs of customers since 2010 with an every growing range of healthcare equipment and the **IKON Mobility** range will allow us to service & supply the needs of the market with the latest battery powered electric mobility products available.

Redgum Brand has been associated with electric mobility devices for nearly 30 years through its parent company FOR-DE Group Pty Ltd, where they have been distributing a leading mobility device brand exclusively during this period. We have a very experienced team that know battery powered mobility and have the resources to back up the **IKON Mobility** range as we increase the number of exciting models in the near future.

With warehouse locations in Perth & Sydney and a dedicated sales & service team **IKON Mobility** by Redgum Brand are ready to assist with your battery powered electric mobility needs.

Visit:



www.ikonmobility.com.au

Contact us:



salesupport@for-de.com.au



08 9248 4180

ASK FOR THE BRAND



Power Mobility

- Enduro Powerlite Pro Front Folding Powered Wheelchair
- Enduro Powerlite Side Folding Electric Wheelchair
- Enduro Powerlite Side Folding Seat Overlays
- Enduro Powerlite Front Folding Electric Wheelchair
- Enduro HD Bariatric Electric Wheelchair
- Enduro Electric Wheelchair Accessories
- STAR Tilt Powerchair
- IKON eMETRO Folding Mobility Scooter
- IKON Move Portable Scooter



ENDURO **POWERLITE PRO**

Carbon Fibre Front Folding Powered Wheelchair




Enduro Powerlite Pro Front Folding Powered Wheelchair

IK2PCXR (Red)



IK2PCXB (Blue)



Remote Control
Included



Turn the page to check out
the Specs and Features >

ENDURO POWERLITE PRO

Features

For the full list of features please refer to the product sheet on the website

- Carbon fibre frame
- Folds up for transport and storage
- Comfortable padded seat
- Remote control included
- Variable speed setting
- Controller can be attached to the **left** or **right** side
- Rear wheel drive
- Controller comes with built in **LED Light**
- Anti-tip wheels
- Fold up** footrest
- Compliant to **ISO7176**



ENDURO POWERLITE PRO

Specifications

For the full list of Specifications please refer to the product sheet on the website

Overall:	970mm (L) x 620mm (W) x 970mm (H)
Overall folded:	630mm (L) x 300mm (W) x 730mm (H)
Seat:	400mm (W) x 425mm (D) x 490mm (H - from floor)
Backrest:	435mm (W) x 525mm (H - from seat)
Speed:	0 - 6km/h
Max. range:	Up to 18km
Product weight:	13.8kg (without battery, cushion and controller)
Max. user weight:	150kg
Motor:	2 x 250W high torque brushless
Battery:	24V / 10Ah Lithium-Ion (supplied)

Additional batteries are available for purchase (sold separately):

Lithium battery replacement for front fold powerlite Pro IK2PCX 10AH - 950100



ENDURO POWERLITE

SIDE FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR




Enduro Powerlite Side Folding Electric Wheelchair

IK2PCS (BLACK)

IK2PCSB (BLUE)

IK2PCSR (RED)

IK2PCSP (PURPLE)



**Remote Control
Included**

Turn the page to check out
the Specs and Features >

ENDURO POWERLITE

SIDE FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR

Features

For the full list of features please refer to the product sheet on the website

- ⦿ **Carbon fibre** frame
- ⦿ **Folds up** for transport and storage
- ⦿ **Comfortable** padded seat
- ⦿ Variable speed setting
- ⦿ Can be used as a **Transit Wheelchair** OR a **Powerchair**
- ⦿ Controller can be attached to the **left** or **right** side
- ⦿ Anti-tip wheels
- ⦿ **Remote control** included
- ⦿ Horn
- ⦿ Compliant to **ISO7176**



ENDURO POWERLITE

SIDE FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR

Specifications

For the full list of Specifications please refer to the product sheet on the website

Overall:	1010mm (L) x 530mm (W) x 895mm (H)
Overall folded:	770mm (L) x 300mm (W) x 740mm (H)
Seat:	490mm (H - from floor)
Between armrests:	460mm
Speed:	0 - 6km/h
Max. range:	Up to 18km
Product weight:	15.8kg
Max. user weight:	140kg
Motor:	2 x 250W high torque brushless
Battery:	24V / 10Ah Lithium-Ion (supplied)

2 battery holders built in allows for a spare battery to be carried (**sold separately**).
Only one battery can operate at a time.

Lithium battery replacement for side fold powerlite IK2PCS 10AH - 950002



Enduro Powerlite Side Folding Electric Wheelchair Seat Overlay Kits

950016 (BLACK OVERLAY KIT)



950017 (BLUE OVERLAY KIT)



950018 (PURPLE OVERLAY KIT)



950019 (RED OVERLAY KIT)



Features

- ⦿ **Extra covers** are now available for purchase to replace on your unit
- ⦿ Attached with velcro strips to allow for **easy removal and attachment**
- ⦿ **Overlay kit includes:**
 - Back cover
 - Seat cover
 - Leg strap

Only fits IK2PCS




ENDURO POWERLITE

FRONT FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR



Enduro Powerlite Front Folding Electric Wheelchair

IK2PCF (BLACK) ■



Remote Control
Included



Turn the page to check out
the Specs and Features >

ENDURO POWERLITE

FRONT FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR

Features

For the full list of features please refer to the product sheet on the website

- Carbon fibre frame
- Folds up for transport and storage
- Comfortable padded seat and high backrest
- Variable speed setting
- Can be used as a **Transit Wheelchair** OR a **Powerchair**
- Storage pocket on backrest
- Anti-tip wheels
- Remote control included
- Controller can be attached to the **left** or the **right** side
- Compliant to **ISO7176**



ENDURO POWERLITE

FRONT FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR

Specifications

For the full list of Specifications please refer to the product sheet on the website

Overall:	920mm (L) x 640mm (W) x 970mm (H)
Overall folded:	720mm (L) x 625mm (W) x 430mm (H)
Seat:	470mm (H - from floor)
Between armrests:	465mm
Speed:	0 - 6km/h
Max. range:	Up to 18km
Product weight:	14.5kg (without battery (1.5kg) and cushion)
Max. user weight:	150kg
Motor:	2 x 250W high torque brushless
Battery:	24V / 10Ah Lithium-Ion (supplied)

Additional batteries are available for purchase (sold separately):

Lithium battery replacement for front fold powerlite IK2PCF 10AH - 950100





ENDURO **POWER HD**

BARIATRIC POWER WHEELCHAIR



Enduro PowerHD Front Folding HD Powered Wheelchair

IK2PCFHD (BLACK) ■



**Remote Control
Included**



*Turn the page to check out
the Specs and Features >*



Features

For the full list of features please refer to the product sheet on the website

- ⦿ **220kg** load capacity
- ⦿ Folds up for **transport** and **storage**
- ⦿ **Comfortable** padded seat
- ⦿ Large joystick knob for **ergonomic preference** due to reduced dexterity
- ⦿ Can be used as a **Transit Wheelchair** OR a **Powerchair**
- ⦿ Anti-tip wheels
- ⦿ **Remote control** included
- ⦿ Front wheel suspension
- ⦿ Controller can be attached to the **left** or the **right** side
- ⦿ **Storage** pocket on backrest



Specifications

For the full list of Specifications please refer to the product sheet on the website

Overall:	1120mm (L) x 720mm (W) x 1020mm (H)
Overall folded:	730mm (L) x 450mm (W) x 830mm (H)
Seat:	500mm (W) x 460mm (D) x 520mm (H - from floor)
Between armrests:	575mm
Speed:	0 - 6km/h
Max. range:	Up to 15km
Product weight:	21.7kg (without battery (2.8kg), cushion & controller))
Max. user weight:	220kg
Motor:	2 x 250W high torque brushless
Battery:	24V / 10Ah Lithium-Ion (supplied)

Additional batteries are available for purchase (sold separately):

Lithium battery replacement for front fold HD powerlite IK2PCF 10AH - 950100





ENDURO ACCESSORIES



Armrest Storage pocket

950014

Features

- ⌚ Space for your everyday items
- ⌚ Easy to fit and remove
- ⌚ Left and right side adaptable

Specifications

Overall: 300mm (W) x 200mm (H)



Bottle Holder

950012

Features

- ⌚ Store your water bottle while on the move
- ⌚ Left and right side adaptable

Specifications

Overall: 118mm (W) x 140mm (H)



Side Storage Pocket for IK2PCS

950013

Features

- ⌚ Space for your everyday items
- ⌚ Easy to fit and remove with two velcro straps
- ⌚ Left and right side adaptable

Specifications

Overall: 225mm (W) x 280mm (H)



Plug In LED Light USB Power Port

X-LUSB

Features

- ⦿ Bright LED light
- ⦿ Dual USB charging ports (A + C)
- ⦿ Plugs into joystick

Attendant Joystick Mount

950023

Features

- ⦿ Mounted to the backrest to allow for attendant control
- ⦿ Fits both IK2PCS & IK2PCF

Specifications

Overall: 155mm (L) x 95mm (W) x 50mm (H)

Coloured Hub Caps

950038 (RED)



950039 (BLUE)



950040 (PURPLE)



Features

- ⦿ Only fits to IK2PCS
- ⦿ Further customise your Enduro Side
- ⦿ Folding Electric Wheelchair



Compact Rear Attendant Controller

950010

Features

- ⦿ Alternative to the remote joystick controller, now in a compact and attachable format
- ⦿ Simple and easy to use
- ⦿ Compatible with IK2PCS, IK2PCF and IK2PCFHD units to allow the attendant to drive from the rear of the unit

Specifications

Overall: 87mm (W) x 60mm (L)



Enduro Swing Away Joystick Mount

950007SW

Features

- ⦿ Allows the joystick controller to swing away from in front of the armrest.
- ⦿ Allows access to end of armrest which makes it easier to get on and off the wheelchair.
- ⦿ Allows wheelchair to get closer to tables.
- ⦿ Left and right versions available.



STAR★ POWERCHAIR



Crash Tested

ISO 7176-19:2022

STAR Tilt Powerchair

IK6PCT (Orange) 



*Turn the page to check out
the Specs and Features >*



Features

For the full list of features please refer to the product sheet on the website

- ⦿ **Ergonomic countoured seat-pad** design to ensure **maximum comfort**
- ⦿ **Pan seat** arrangement for various seating systems for **maximum adaptability**
- ⦿ **Advanced LINX programmable joystick controller system** provides the user with both **independence and safety**
- ⦿ **Adjustable width** seat & backrest
- ⦿ **Power recline** function
- ⦿ **Power tilt** seating
- ⦿ **Adjustable** footrests, armrests and headrest
- ⦿ **Crash Tested** and can be used as a seat in a transport vehicle



Specifications

For the full list of Specifications please refer to the product sheet on the website

Overall:	1080mm (L) x 580mm (W) x 1215mm (H)
Overall with legrest:	1155 - 1820mm (L)
Seat:	430 - 480mm (D) x 430 - 480mm (W)
Seat tilt range:	0° - 25°
Backrest angle range:	90° - 140°
Speed:	0 - 10km/h
Max. range:	Up to 30km
Product weight:	106.5kg
Max. user weight:	136kg



ISO 7176-19:2022



emetro
portable mobility scooter



IKON eMETRO Portable Mobility Scooter

IKSSER (RED)



IKSSEY (YELLOW)



IKSSEB (BLUE)



Turn the page to check out
the Specs and Features >



Features

For the full list of features please refer to the product sheet on the website

- ⓘ **Automatically slows** when turning
- ⓘ **Folds up for transport and storage**
- ⓘ **Ergonomic** padded seat
- ⓘ **Pneumatic** tyres
- ⓘ **Ergonomic** hand grips
- ⓘ **Anti-tip** wheels
- ⓘ **Removable** carry bag
- ⓘ **Flip down** armrest for side transfer



**Aeroplane
Travel Ready**



Specifications

For the full list of Specifications please refer to the product sheet on the website

Overall:	970mm (L) x 540mm (W) x 900mm (H)
Overall folded:	970mm (L) x 540mm (W) x 310mm (H)
Seat:	400mm (W) x 325mm (D) x 405mm (H - from footrest)
Backrest height:	310mm
Speed:	0 - 8km/h
Max. range:	Up to 12km
Product weight:	18.6kg (excl. battery *(2kg))
Max. user weight:	120kg

Additional batteries are available for purchase (sold separately):

10AH Lithium battery replacement for eMETRO IK5SE - 960009

12AH Lithium battery replacement for eMETRO IK5SE - 960010



MOVE

MOBILITY SCOOTER



IKON Move Portable Mobility Scooter

IK4SEP (METALLIC PURPLE)



IK4SEB (METALLIC BLUE)



IK4SEO (GOLDEN ORANGE)



Turn the page to check out
the Specs and Features >

MOVE

PORTABLE MOBILITY SCOOTER

Features

For the full list of features please refer to the product sheet on the website

- ⌚ **Pulls apart** for transport and storage
- ⌚ **Comfortable** contoured seat
- ⌚ **Full** suspension
- ⌚ **Padded** loop handle bar
- ⌚ **Angle adjustable** tiller
- ⌚ **Removable** shopping basket
- ⌚ **Bright** LED headlight
- ⌚ **Rear** wheel drive



MOVE

PORTABLE MOBILITY SCOOTER

Specifications

For the full list of Specifications please refer to the product sheet on the website

Overall length:	1080mm
Overall width:	590mm
Overall height:	880mm
Max. range:	Up to 18km
Speed:	0 - 6.4km/h
Max. range:	Up to 18km
Batteries:	2 x 12V / 20Ah SLA (supplied)
Product weight:	58kg
Max. user weight:	130kg





Bathroom aids. Mobility Aids. Ramps.

We do it all. Just *ask for the brand.*

Check out our Redgum Brand Range at

www.redgumbrandmobility.com.au



Scooters. Powerchairs.

We do it all. Just *ask for the brand.*

Check out our Shoprider Range at

www.shoprider.com.au



REDGUM Brand: Various Assisted Living and Lifestyle Products

SHOPRIDER Mobility: Electric Mobility Scooters and Powerchairs

IKON Scooters: Electric Mobility Scooters and Powerchairs

Scooter Xtras: Mobility Scooter and Powerchair Accessories

Baricare: Various Bariatric Assisted Living and Lifestyle Products

IKON Mobility owns the copyright to this product catalogue. Any information contained in this catalogue shall not be photocopied, reproduced or translated into any other languages without prior written consent of **IKON Mobility**.

IKON Mobility strives to provide the highest quality products possible that are subject to change without notice.

© Copyright **IKON Mobility**

IKON Mobility:

Registered NDIS Provider

TGA (Therapeutic Goods Administration) Registered

Proudly 100% Australian Owned and Operated

ATSA Member

Perth / Sydney Offices + Warehousing

Distributors Australia wide

In-house R+D team - Product Development

In-house servicing of all products

Disclaimer:

The lifestyle images in this catalogue have been made using AI. This may cause the product to differ slightly to what it actually is.



ASK FOR THE BRAND

