



ENDURO POWERLITE

FRONT FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR

IK2PCF

CARBON FIBRE



14.5kg
NET WEIGHT

BLACK
IK2PCF



ISO7176

ISO9999 : 12 23 03

Electrically powered wheelchairs. Products
for personal mobility and transportation

FEATURES

- Comfortable padded seat and backrest
- Storage pocket on backrest
- Onboard controller with LCD display and joystick controls
- Remote control included
- Battery power indicator
- Onboard controller can be attached on left or right hand side
- Variable speed setting
- Horn
- Rear wheel drive
- Front wheel suspension
- Electromagnetic brake
- Airline safe battery
- Rechargeable Lithium-ion battery
- Anti-tip wheels

P: **08 9248 4180**
E: **orders@for-de.com.au**



Scan the QR code for more
information on this product

Please turn over for
further information

ASK FOR THE BRAND

ikonmobility.com.au



ENDURO POWERLITE

FRONT FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR

IK2PCF

CARBON FIBRE



CARBON FIBRE
Frame

CONVENIENT
Storage Bag



**AIRLINE
TRAVEL
READY**



FOLD UP
For Transport
and Storage

FOLD DOWN
With one Movement

FOLDING
Foot Plate



Drive/Manual
Adjustment

**MOTORISED OR
FREEWHEEL**

CONTROLLER
Onboard + Remote



Anti Tip
Wheels

ASK FOR THE BRAND

ikonmobility.com.au



ENDURO POWERLITE

FRONT FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR

IK2PCF


SPECIFICATIONS

All dimensions are nominal and may differ +/- 10mm




Overall length (L)

920mm




Seat depth (L)

420mm




Push handle height (H)

970mm




Armrest length (L)

390mm




Between armrests (W)

460mm




Safe climb angle

13°




Overall height (H)

970mm




Seat height (H - from floor)

510mm



Backrest height (H - from seat)

520mm



Armrest height (H - from seat)

200mm




Product weight

**14.5kg net weight
16.2kg gross weight**




Speed

6km/h




Overall width (W)

640mm




Backrest width (W)

450mm



Seat width (W)

420mm




Armrest width (W)

50mm



Max. user weight

150kg (on flat ground)



Max. range

Up to 18km

ASK FOR THE BRAND

ikonmobility.com.au



ENDURO POWERLITE

FRONT FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR


IK2PCF

SPECIFICATIONS


All dimensions are nominal and may differ +/- 10mm




Front tyres (flat free)
165mm (6.5")




Tipping wheel distance
(from floor) **30mm**




Footplate to seat (D)
440mm



Motor
**2 x 250W high torque
brushless motors**




Rear tyres (flat free)
Rear: 280mm (11")




Footplate length
300mm




Footplate angle
90°



Battery
24V / 10Ah Lithium-Ion




Push handles width (W)
440mm




Footplate width
145mm




Charger
2 amp Lithium-Ion



Overall length (L)
970mm



Overall height (H)
430mm



Overall width (W)
340mm

FOLDED

ASK FOR THE BRAND

ikonmobility.com.au