

ENDURO POWER HD

BARIATRIC CARBON FIBRE POWER WHEELCHAIR

Please turn over for more information



ASK FOR THE BRAND

ikonmobility.com.au

ENDURO POWER HD

FRONT FOLDING POWERED WHEELCHAIR

220KG Max.
user weight

500mm wide seat

Folds for storage
& transport

IK2PCFHD

Please contact **IKON Mobility** if you require further information on this product

P: **08 9248 4180**
E: **orders@for-de.com.au**

FEATURES

- 220kg load capacity
- Carbon fibre frame
- Folds up for transport and storage
- Comfortable padded seat and backrest
- Lift up armrests
- Onboard controller with LCD display and fingertip controls
- Remote control included
- Battery power indicator
- Controller can be attached left or right side
- Variable speed setting
- Rear wheel drive
- Front wheel suspension
- Electromagnetic brake
- Rechargeable Lithium-ion battery
- Aeroplane safe battery
- Fold up footrests
- Anti-tip wheels
- Lap belt
- Horn

SPECIFICATIONS

- **Overall:** 1120mm (L) x 720mm (W) x 1020mm (H)
- **Overall folded:** 730mm (L) x 450mm (W) x 830mm (H)
- **Seat:** 500mm (W) x 460mm (D) x 520mm (H from floor)
- **Backrest:** 550mm (W) x 475mm (H from seat)
- **Push handle height:** 1020mm
- **Turning circle:** 860mm
- **Speed:** 0 - 6km/h
- **Maximum range:** Up to 15km*
- **Safe climb angle:** 0 - 10°
- **Maximum climb angle:** 13°
- **Motor:** 2 x 250W high torque brushless motors
- **Battery:** 24V / 10Ah Lithium-Ion (supplied)
- **Charger:** 2.0 Amp Lithium-Ion (Supplied)
- **Tyres (Flat free):** Front: Ø200mm / Rear: Ø325mm
- **Product weight:** 21.7kg without battery, cushion and controller
- **Battery weight:** 2.8kg
- **Maximum user weight:** 220kg (on flat ground)

ASK FOR THE BRAND

ikonmobility.com.au