

ENDURO POVERLITE FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR

CARBON FIBRE



FEATURES

- Carbon fibre frame
- Folds up for transport and storage
- Comfortable padded seat and backrest
- Fold down backrest with storage pocket
- Lift up armrests
- Onboard controller with LCD display and fingertip controls
- Battery power indicator
- Controller can be attached left or right side
- Variable speed setting
- Rear wheel drive
- Front wheel suspension
- Electromagnetic brake
- Aeroplane safe battery
- Rechargeable Lithium-ion battery
- Anti-tip wheels
- Calf strap
- Lap belt
- Horn





ENDURO POVERLITE IKZPCS FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR

CARBON FIBRE



SPECIFICATIONS

- Overall (L to unfolded footrest):1045mm(L) x 555mm(W) x 885mm(H)
- Folded: 820mm(L) x 275mm(w) x 770mm(H)
- Seat (H from floor):420mm(L) x 460mm(W) x 490mm(H)
- Backrest (H from seat): 440mm(W) x 420mm (H)
- Push Handle Height: 880mm
- Turning Circle: 1080mm
- Speed: **6km/hr**
- Max Range: **Up to 18km**¹
- Safe Climb Angle: **6 10**°
- Max Climb Angle: 10°
- Motor: **2 x 250W high torque brushless motors**
- Battery: 24V / 10Ah Lithium-Ion
- Charger: **2 Amp Lithium-Ion**
- Tyres (flat free):Front: 165mm (6.5") / Rear: 280mm (11")
- Product Weight:13.2kg net weight,15.8kg gross weight
- Max User Weight: 140kg (on flat ground)

¹Dependant on occupant and conditions.

ASK FOR THE BRAND

ikonmobility.com.au