



ENDURO POWERLITE^{IK2PCS}

SIDE FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR

CARBON FIBRE



13.2kg
NET WEIGHT

BLACK
IK2PCS



BLUE
IK2PCSB



PURPLE
IK2PCSP



RED
IK2PCSR

FEATURES

- Carbon fibre frame
- Folds up for transport and storage
- Comfortable padded seat and backrest
- Fold down backrest with storage pocket
- Lift up armrests
- Onboard controller with LCD display and joystick controls
- Battery power indicator
- Controller can be attached left or right side
- Variable speed setting
- Rear wheel drive
- Front wheel suspension
- Electromagnetic brake
- Aeroplane safe battery
- Rechargeable Lithium-ion battery
- Anti-tip wheels
- Calf strap
- Lap belt
- Horn
- 3 Different Colours



**AIRLINE
TRAVEL
READY**

Please turn over for
product specifications

ASK FOR THE BRAND

ikonmobility.com.au



ENDURO POWERLITE

SIDE FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR

CARBON FIBRE



**MOTORISED OR
FREEWHEEL**

**ANTI-TIP
WHEELS**

**CONVENIENT
STORAGE BAG**

**FLIP UP
ARMRESTS**

**FOLDS UP
FOR STORAGE**



**REMOTE
CONTROL
AVAILABLE**



SPECIFICATIONS

- Overall (L to unfolded footrest):
1045mm(L) x 555mm(W) x 885mm(H)
- Folded:
820mm(L) x 275mm(w) x 770mm(H)
- Seat (H from floor):
420mm(L) x 460mm(W) x 490mm(H)
- Backrest (H from seat):
440mm(W) x 420mm (H)
- Push Handle Height: **880mm**
- Turning Circle: **1080mm**
- Speed: **6km/hr**
- Max Range: **Up to 18km¹**
- Safe Climb Angle: **6 - 10°**
- Max Climb Angle: **10°**
- Motor: **2 x 250W high torque brushless motors**
- Battery: **24V / 10Ah Lithium-Ion**
- Charger: **2 Amp Lithium-Ion**
- Tyres (flat free):
Front: 165mm (6.5") / Rear: 280mm (11")
- Product Weight:
**13.2kg net weight,
15.8kg gross weight**
- Max User Weight: **140kg (on flat ground)**

P: **08 9248 4180**
E: **orders@for-de.com.au**



For more information
about this product

ASK FOR THE BRAND

ikonmobility.com.au