

## SIDE FOLDING



Hammock cushion • Handles for pushing • Double plate footrest



## FRONT FOLDING



Solid cushion • Crossbar for pushing • Single plate footrest



## SPECIFICATIONS

(Specs may differ  
+/- 10mm)

### SIDE

### FRONT

#### Overall folded:

**L:** 820mm  
**W:** 275mm  
**H:** 770mm

**L:** 630mm  
**W:** 340mm  
**H:** 720mm

#### Seat:

(H from floor to  
base of cushion)

**D:** 420mm  
**W:** 460mm  
**H:** 470mm

**D:** 400mm  
**W:** 400mm  
**H:** 460mm

#### Backrest:

(H from seat)

**W:** 440mm  
**H:** 420mm

**W:** 440mm  
**H:** 500mm

#### Push handle height:

880mm

970mm

#### Turning circle:

1080mm

750mm

#### Safe climb angle:

0 - 10°

0 - 13°

#### Max. climb angle:

10°

13°

#### Product weight: (without battery)

13.2kg

14.5kg

#### Max. user weight:

140kg

150kg