

ENDURO POWERLITE PRO

FRONT FOLDING POWERED WHEELCHAIR

Please turn over for more information



ASK FOR THE BRAND

ikonmobility.com.au



IK2PCX

FEATURES

- Carbon fibre frame
- Folds up for transport and storage
- Comfortable padded seat
- Mesh backrest
- Lift up armrests
- Onboard controller with LCD display and fingertip controls
- Remote control included
- Battery power indicator
- Controller can be attached left or right side
- Variable speed setting
- Rear wheel drive
- Front wheel suspension
- Electromagnetic brake
- Rechargeable Lithium-ion battery
- Aeroplane safe battery
- Fold up footrest
- Anti-tip wheels
- Lap belt
- Horn

SPECIFICATIONS

- **Overall:** 970mm (L) x 620mm (W) x 970mm (H)
- **Overall Folded:** 630mm (L) x 300mm (W) x 730mm (H)
- **Seat:** 400mm (W) x 425mm (D) x 490mm (H from floor)
- **Backrest:** 435mm (W) x 525mm (H from seat)
- **Turning Circle:** 750mm
- **Speed:** 0 - 6km/h
- **Maximum range:** Up to 18km*
- **Safe climb angle:** 0 - 10°
- **Maximum climb angle:** 10°
- **Motor:** 2 x 250W high torque brushless motors
- **Battery:** 24V / 10Ah Lithium-Ion (supplied)
- **Charger:** 2.0 Amp Lithium-Ion (Supplied)
- **Tyres (Flat free):** Front: Ø160mm / Rear: Ø210mm
- **Product Weight:** 13.8kg (without battery, cushion and controller)
- **Battery weight:** 1.8kg
- **Maximum User Weight:** 130kg (on flat ground)

* Maximum range is based on an ambient temperature of 20°C, 75kg user, brand new fully charged battery and constant driving speed of 5km/h on flat level ground.

Please contact **IKON Mobility** if you require further information on this product:

P: **08 9248 4180**
E: **orders@for-de.com.au**