



Part of the **FOR-DE Group Pty Ltd**

3 NEW COLOURS

GOLDEN ORANGE
IK4SEO

METALLIC BLUE
IK4SEB

METALLIC PURPLE
IK4SEP



MOVE

PORTABLE MOBILITY SCOOTER

ASK FOR THE BRAND

ikonmobility.com.au

IMPORTANT NOTE: It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

IKON MOBILITY strives to provide the highest quality products possible that are subject to change without notice.



Part of the **FOR-DE Group Pty Ltd**

NEW COLOURS

PULLS APART for transport and storage

**GOLDEN
ORANGE**
IK4SEO

**METALLIC
BLUE**
IK4SEB

**METALLIC
PURPLE**
IK4SEP

MOVE

PORTABLE MOBILITY SCOOTER

FEATURES

- Portable scooter
- Pulls apart for transport and storage
- Comfortable mid back contoured seat
- Full improved suspension
- Padded loop handlebar
- Angle adjustable tiller
- Width adjustable flip up padded armrests
- Removable shopping basket
- LED Headlight
- Horn
- Fingertip operated controls
- Forward / Reverse drive settings
- Variable speed setting
- Battery power indicator
- Rear wheel drive
- Flat free tyres
- Electromagnetic brake
- Front nudge bar
- Anti tip wheels at rear

WARRANTY

- Structural Frame: **5 years**
- Drive Train Assembly: **18 months**
- Other Parts: **12 months**
- 20 Amp SLA Batteries: **18 months**
- 4 Amp Charger: **18 months**

Please turn over for product specifications



PORTABLE SCOOTER

P: **08 9248 4180**

E: **orders@for-de.com.au**

ASK FOR THE BRAND

ikonmobility.com.au

IMPORTANT NOTE: It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

IKON MOBILITY strives to provide the highest quality products possible that are subject to change without notice.



OVERALL LENGTH

1080mm



OVERALL WIDTH

590mm



OVERALL HEIGHT

880mm



FRONT ASSEMBLY

L: 940mm
W: 450mm
H: 410mm
Weight: 19kg



REAR WHEEL ASSEMBLY

L: 320mm
W: 49mm
H: 250mm
Weight: 13kg



BATTERY PACK

W: 330mm
D: 230mm
H: 230mm
Weight: 14kg




SEAT (FOLDED DOWN)

W: 460mm
D: 610mm
H: 530mm
Weight: 11kg




BASKET

W: 270mm
D: 210mm
H: 240mm
Weight: 1kg



SEAT WIDTH (INT)

450mm



SEAT DEPTH

440mm



SEAT HEIGHTS

330mm
335mm
380mm



ARMREST LENGTH

320mm



ARMREST WIDTH

80mm



ARMREST HEIGHT

200mm




TILLER WIDTH

400mm



TILLER HEIGHT

760mm



TURNING CIRCLE

1400mm



SPEED

6.4km / h



MAXIMUM RANGE

Up to 18km



MAXIMUM CLIMB ANGLE

10°



MOTOR

270w high torque motor



BATTERY PACK

2x 12V / 20 Ah
(supplied)



BATTERY CHARGER

4 Amp
(supplied)



TYRES (FLAT FREE)

Ø200mm



PRODUCT WEIGHT

58kg



MAXIMUM USER WEIGHT

130kg
(on flat ground)

Please contact IKON Mobility if you would like any further information on the Move Portable Mobility Scooter.

ASK FOR THE BRAND

ikonmobility.com.au