



PULLS APART
for transport and storage

Part of the **FOR-DE Group Pty Ltd**

NEW COLOUR!



**PORTABLE
SCOOTER**



METALLIC BLUE
IK4SEB

METALLIC ORANGE
IK4SEO

METALLIC GREEN
IK4SEG

MOVE

PORTABLE MOBILITY SCOOTER

FEATURES

- Portable scooter
- Pulls apart for transport and storage
- Comfortable mid back contoured seat
- Full suspension
- Padded loop handlebar
- Angle adjustable tiller
- Width adjustable flip up padded armrests
- Removable shopping basket
- LED Headlight
- Horn
- Fingertip operated controls
- Forward / Reverse drive settings
- Variable speed setting
- Battery power indicator
- Rear wheel drive
- Flat free tyres
- Electromagnetic brake
- Front nudge bar
- Anti tip wheels at rear

WARRANTY

- Structural Frame: **5 years**
- Drive Train Assembly: **18 months**
- Other Parts: **12 months**
- 20 Amp Hybrid Gel Batteries: **18 months**
- 2 Amp Charger: **18 months**

Flip over for Specifications

P: **08 9248 4180**
E: **orders@for-de.com.au**

ASK FOR THE BRAND ikonmobility.com.au

IMPORTANT NOTE: It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.




















IKON MOBILITY strives to provide the highest quality products possible that are subject to change without notice.



Please contact IKON Mobility if you would like any further information on the Move Portable Mobility Scooter.

MOVE

PORTABLE MOBILITY SCOOTER

								
OVERALL LENGTH	OVERALL WIDTH	OVERALL HEIGHT	FRONT ASSEMBLY	REAR WHEEL ASSEMBLY	BATTERY PACK	SEAT (FOLDED DOWN)	BASKET	SEAT WIDTH (INT)
1080mm	590mm	880mm	L: 940mm W: 450mm H: 410mm Weight: 19kg	L: 320mm W: 49mm H: 250mm Weight: 13kg	W: 330mm D: 230mm H: 230mm Weight: 16kg	W: 460mm D: 610mm H: 530mm Weight: 11kg	W: 270mm D: 210mm H: 240mm Weight: 1kg	450mm
								
SEAT DEPTH	SEAT HEIGHTS	ARMREST LENGTH	ARMREST WIDTH	ARMREST HEIGHT	TILLER WIDTH	TILLER HEIGHT	TURNING CIRCLE	SPEED
440mm	330mm 335mm 380mm	320mm	80mm	200mm	400mm	760mm	1400mm	6.4km / h
								
MAXIMUM RANGE	SAFE CLIMB ANGLE	MAXIMUM CLIMB ANGLE	MOTOR	BATTERY PACK	BATTERY CHARGER	TYRES (FLAT FREE)	PRODUCT WEIGHT	MAXIMUM USER WEIGHT
Up to 18km	6°	10°	270w high torque motor	2x 24V / 20 Ah (supplied)	2 Amp (supplied)	Ø200mm	60kg	130kg (on flat ground)

ASK FOR THE BRAND

ikonmobility.com.au

IMPORTANT NOTE: It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

IKON MOBILITY strives to provide the highest quality products possible that are subject to change without notice.