



Part of the **FOR-DE Group Pty Ltd**



**IMPORTANT NOTE:**

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.



**PULLS APART**  
for transport and storage



METALLIC BLUE  
IK4SEB

METALLIC ORANGE  
IK4SEO

# MOVE

PORTABLE MOBILITY SCOOTER

## FEATURES

- Portable scooter
- Pulls apart for transport and storage
- Comfortable mid back contoured seat
- Full suspension
- Padded loop handlebar
- Angle adjustable tiller
- Width adjustable flip up padded armrests
- Removable shopping basket
- LED Headlight
- Horn
- Fingertip operated controls
- Forward / Reverse drive settings
- Variable speed setting
- Battery power indicator
- Rear wheel drive
- Flat free tyres
- Electromagnetic brake
- Front nudge bar
- Anti tip wheels at rear

## WARRANTY

- Structural Frame: **5 years**
- Drive Train Assembly: **18 months**
- Other Parts: **12 months**
- XT 20 Amp Hybrid Gel Batteries: **18 months**
- 2 Amp Charger: **18 months**

**Flip over for Specifications**

P: **08 9248 4180**

E: **orders@for-de.com.au**

ASK FOR THE BRAND

ikonmobility.com.au

IKON Mobility strives to provide the highest quality products possible that are subject to change without notice.



Please contact IKON Mobility if you would like any further information on the Move Portable Mobility Scooter.

# MOVE

PORTABLE MOBILITY SCOOTER

<b>OVERALL LENGTH</b>	<b>OVERALL WIDTH</b>	<b>OVERALL HEIGHT</b>	<b>FRONT ASSEMBLY</b>	<b>REAR WHEEL ASSEMBLY</b>	<b>BATTERY PACK</b>	<b>SEAT (FOLDED DOWN)</b>	<b>BASKET</b>	<b>SEAT WIDTH (INT)</b>
1080mm	590mm	880mm	L: 940mm W: 450mm H: 410mm Weight: 19kg	L: 320mm W: 49mm H: 250mm Weight: 13kg	W: 330mm D: 230mm H: 230mm Weight: 16kg	W: 460mm D: 610mm H: 530mm Weight: 11kg	W: 270mm D: 210mm H: 240mm Weight: 1kg	450mm
<b>SEAT DEPTH</b>	<b>SEAT HEIGHTS</b>	<b>ARMREST LENGTH</b>	<b>ARMREST WIDTH</b>	<b>ARMREST HEIGHT</b>	<b>TILLER WIDTH</b>	<b>TILLER HEIGHT</b>	<b>TURNING CIRCLE</b>	<b>SPEED</b>
440mm	330mm 335mm 380mm	320mm	80mm	200mm	400mm	760mm	1400mm	6.4km / h
<b>MAXIMUM RANGE</b>	<b>SAFE CLIMB ANGLE</b>	<b>MAXIMUM CLIMB ANGLE</b>	<b>MOTOR</b>	<b>BATTERY PACK</b>	<b>BATTERY CHARGER</b>	<b>TYRES (FLAT FREE)</b>	<b>PRODUCT WEIGHT</b>	<b>MAXIMUM USER WEIGHT</b>
Up to 18km	6°	10°	270w high torque motor	2x 24V / 20 Ah (supplied)	2 Amp (supplied)	Ø200mm	60kg	130kg (on flat ground)

ASK FOR THE BRAND

ikonmobility.com.au