



IKSSEY, IKSSER

eMetro

PORTABLE MOBILITY SCOOTER

Please turn over for more information



ASK FOR THE BRAND

ikonmobility.com.au

IMPORTANT NOTE: It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

IKON strives to provide the highest quality products possible that are subject to change without notice.

Removeable
Storage
Basket



Folds for
Transport &
Storage

IK5SEY IK5SER



Please contact **REDGUM Brand** if you
require further information on this product

P: **08 9248 4180**
E: **orders@for-de.com.au**

FEATURES

- Folds up for transport and storage
- Comfortable padded seat
- Under seat storage
- Removable carry bag
- Mesh backrest
- Width adjustable drop down armrests
- Angle adjustable tiller
- Onboard controller with LCD display and fingertip controls
- Battery power indicator
- Variable speed setting
- Automatically slows when turning
- Ergonomic hand grips
- Rear wheel drive
- Front wheel suspension
- Electromagnetic brake
- Rechargeable Lithium-ion battery
- Aeroplane safe battery
- Anti-tip wheels
- Horn

SPECIFICATIONS

- **Overall:** 970mm (L) x 540mm (W) x 900mm (H)
- **Overall folded:** 970mm (L) x 540mm (W) x 310mm (H)
- **Seat:** 400mm (W) x 325mm (D) x 405mm (H from footrest)
- **Backrest:** 400mm (W) x 310mm (H from seat)
- **Handle height:** 865mm (Angle adjustable)
- **Turning circle:** 1100mm
- **Speed:** 0 - 8km/h
- **Maximum range:** Up to 12km*
- **Safe climb angle range:** 0 - 9°
- **Maximum climb angle:** 9°
- **Motor:** 24v 250w Brushless motor
- **Battery:** 10.4Ah Lithium-Ion (supplied)
- **Charger:** 2Amp Lithium-Ion (supplied)
- **Tyres (Pneumatic):** Ø230mm
- **Product Weight:** 18.6kg without battery 2kg battery
- **Maximum User Weight:** 120kg (on flat ground)

ASK FOR THE BRAND

ikonmobility.com.au