



STAR TILT POWERCHAIR

IK6PCT

STAR★
POWERCHAIR

COMING SOON



FEATURES

- Ergonomic seat design to ensure maximum comfort
- Pan seat arrangement for various seating systems
- Advanced technology provides the user with both independence and safety
- Adjustable width seat & backrest
- Recline function
- Power tilt seating
- High ground clearance with 14" drive wheel
- Certified safety belt
- Adjustable footrests, armrests and headrest
- Switch between driving/seat tilt mode easily
- Front & rear lights
- Rear suspension
- Anti-tip wheels
- Drive & freewheel mode
- Elevating swing away legrests
- Fold up footplates
- Large front wheels

 **CRASH TESTED**
ISO 7176-19



Scan the **QR** code to visit the
IKON mobility website!

P: **08 9248 4180**
E: **orders@for-de.com.au**

ASK FOR THE BRAND

ikonmobility.com.au

Measurements are nominal and may differ +/- 10mm

SPECIFICATIONS

Overall Length: **1080mm**

Overall Width: **580mm**

Overall Height: **1215mm**

Overall (L) with legrest: **1155 - 1820mm**

Folded length: **850mm**

Folded width: **600mm**

Folded height: **970mm**

Backrest height: **520mm**

Backrest width: **430 - 480mm**

Backrest angle range: **90° - 144°**

Footrest to seat: **450 - 575mm**

Leg to seat angle: **110°**

Floor to seat pan: **453mm**

Floor to std. seat cushion: **510mm**

Armrest to seat distance: **170 - 290mm**

Seat depth: **430 - 480mm**

Seat width: **430 - 480mm**

Seat tilt range: **0° - 22.5°**

Motor: **2 x 340 W**

Battery: **24V52A**

Battery charger: **6 Amp**

Rear wheels diameter: **14 x 3 "**

Front wheels diameter: **9 x 3 "**

Max. user weight: **136kg**

Product weight: **106.5kg**

Battery weight: **26kg**

Max. speed: **10km/h**

Min braking distance at max speed: **1700mm**

Obstacle clearance: **50mm**

Safe climbing angle: **10°**

Turning Radius: **900mm**

Range: **30km**



CRASH TESTED
ISO 7176:19



P: **08 9248 4180**
E: **orders@for-de.com.au**

ASK FOR THE BRAND

ikonmobility.com.au