



METALLIC SILVER IK8PRS

PRESTO MOBILITY SCOOTER

Part of the FOR-DE Group Pty Ltd



EUROPEAN DESIGN



FEATURES

- Latest European design features
- High back contoured seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- LED Headlight
- Rear lights
- LCD display
- Indicators / Hazard lights
- Large battery capacity
- Battery power indicator
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width armrests
- Removable shopping basket

WARRANTY

Structural Frame: **5 years**

Drive Train Assembly: **18 months**

Other Parts: **12 months**

XT 75 Amp Hybrid Gel Batteries: **18 months**

Victron Premium Charger: **5 years**

Flip over for Specifications

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

ASK FOR THE BRAND

ikonmobility.com.au



SPECIFICATIONS


PRESTO MOBILITY SCOOTER



OVERALL LENGTH
1340mm




OVERALL HEIGHT
1300mm



OVERALL WIDTH
740mm




SEAT WIDTH
510mm



SEAT DEPTH
480mm



SEAT HEIGHT
450 - 525mm
(4 adjustments)



ARMREST LENGTH
380mm



ARMREST WIDTH
60mm




**ARMREST HEIGHT
(FROM SEAT)**
190mm



TILLER WIDTH
470mm



TILLER HEIGHT
930mm



TURNING CIRCLE
1600mm



SPEED
10km/h



MAXIMUM RANGE
Up to 44km




SAFE CLIMB ANGLE
10°




MAXIMUM CLIMB ANGLE
15°



MOTOR
2.0Hp high torque motor
(4 Pole)



BATTERY
2 x 75Amp minimum
(supplied)



BATTERY CHARGER
8 Amp



TYRES (PNEUMATIC)
Ø356mm - Front + Rear



PRODUCT WEIGHT
136kg



MAXIMUM USER WEIGHT
160kg
(on flat ground)

Please contact IKON Mobility if you require more information

P: **08 9248 4180**

E: **orders@for-de.com.au**

ASK FOR THE BRAND

ikonmobility.com.au